



BIOTRONIK


CARDIOMESSENGER SET UP AND USE

Setting up BIOTRONIK Home Monitoring is simple and, once it is set up, the system is fully automatic. As a patient, you are not required to take any action for successful daily transmissions to your care team.

You will receive a CardioMessenger device either at the time of your procedure or shortly thereafter in the mail. It is important to follow the instructions given to you at the time of your procedure or refer to the Quick Start Guide included with your CardioMessenger. Together with your physician, you can decide at what time of the day data transmission will take place. Typically, these transmissions are scheduled while you sleep.

After returning home after your procedure, please set up the CardioMessenger as soon as possible, preferably the same day. Setup is easy, with just three simple steps.

- Step 1: Remove contents from CardioMessenger box.**
- Step 2: Position the CardioMessenger near your bed.**
- Step 3: Plug the CardioMessenger into a nearby power outlet.**



After it is plugged in, the CardioMessenger will turn on automatically and performs a self-test. The CardioMessenger is ready for use once the self-test is completed and the following icons are displayed:



If this is not the case, please refer to the CardioMessenger Smart Patient Manual. Both the CardioMessenger Smart Quick Start Guide and the CardioMessenger Smart Patient Manual contain important information concerning the setup and use of your CardioMessenger. Please review each thoroughly.



If you need assistance please contact HeartHQ's Remote Monitoring Team on 07 5414 1100.